R

Quantaurus-Tau



Quantaurus-Tau is a compact system for measuring fluorescence lifetimes in the sub nanosecond to millisecond range. Operation is simple, just set the sample into the sample chamber, and enter a few conditions on the measurement software to measure the fluorescence lifetime and PL spectrum in a short time with high precision. In a typical measurement, analysis results are obtained in a mere 60 seconds.

Fluorescence Lifetime

Fluorescence lifetime measurement

TCU

Measuring an excited-state relaxation process

The fluorescence spectrum obtained from an organic material or fluorescent probe is a vital parameter for controlling and evaluating the material functions and characteristics such as the peak wavelength and fluorescence intensity. However, a fluorescence spectrum usually shows time-integrated information, and so when the material contains multiple substances and reactive elements, their fluorescence spectrum can only be acquired as integrated information. An effective approach in such cases is to observe the light emission dynamics by making use of the time axis parameter. This is generally called fluorescence lifetime measurement, in which the time required for the substance excited by the pulsed light to return to its ground state is measured in the sub-nanosecond to millisecond region. This measurement allows obtaining more information such as multiple different fluorescence lifetimes even at the same wavelength and the percentage in which they are present within the material, etc.

Features

- High sensitivity measurement by photon counting method
- Time resolution better than 100 ps (by deconvolution)

11

- Cooling function for solution sample (-196 °C) (option)
- Phosphorescence measurement (option)
- Time-resolved measurement of fluorescence anisotropy (option)
- Fluorescence spectrum measurement
- Space-saving, compact design

Easy and quick measurements

Emission Lifetime can be gotten easily and quickly only by putting the sample into sample box and setting the 4 measurement conditions.

7 excitation wavelength

280 nm, 340 nm, 365 nm, 405 nm, 470 nm, 590 nm, and 630 nm.

Analyzing diffrent sample forms

Thin-film, solid, solutions and powder.

2 selection of detector





Fluorescence lifetime measurement is applicable to varied applications. Typical applications include electron movement and energy transfer reactions within or between organic metallic complex molecules, as well as fluorescence and phosphorescence lifetime measurement of materials essential for developing organic EL devices, FRET (fluorescence resonance energy transfer) in fluorescent proteins, and pass/fail testing of compound semiconductors for solar cells and LED, etc.

NIR

C11367-12/-15

1100



Measurement procedure

Software designed taking account of the measurement procedure ensures easy and quick measurements.



Analysis functions

Quantaurus-Tau includes a variety of measurement and analysis functions such as simultaneous multi-component fluorescence lifetime measurement and multi-sample data comparison.

Multi-component fluorescence lifetime analysis and comparison



Multi-sample fluorescence lifetime analysis and comparison



Multi-sample PL spectrum analysis and comparison



Multi-component analysis of up to five components

In fluorescence lifetime measurement, a phenomenon often occurs where the data is observed as the sum of the attenuation curves of multicomponent fluorescence lifetimes. Quantaurus-Tau easily calculates the fluorescence lifetime data and component ratio of each element by using the dedicated software.

Highly accurate analysis by deconvolution

Deconvolution processing enables fluorescence lifetime analysis with high accuracy. When analyzing longer lifetime components such as phosphorescence, the "Tail Fit" function can be used instead of deconvolution processing.

Real-time display of time profiles and spectrum

Time profiles or spectrum are displayed on the monitor screen in real-time. This is a useful function for selecting the time scale during measurement or determining the analysis data range.

Multiple data analysis on the same screen

Calculated fluorescence lifetime values are also displayed on the same screen for easy comparison analysis.

Comparisons under the same fitting

To make comparison analysis under the same conditions, Quantaurus-Tau subjects the multiple samples to specific fitting ranges, IRF (instrument response functions), and parameter settings.

Graph editing with a graph setup feature

This allows you to change the range of each axis as needed on the comparison analysis screen so that the data can be edited to match your purpose. This feature also allows powerful normalizing whenever needed.

Acquired data can be easily stored as text data

The acquired data can be stored into the graph analysis software as text data by simple copy-and-paste operations.

Time-resolved spectrum display

Allows time-resolved spectrum display the greatest feature offered by streak camera systems.

Spectrum and fluorescence decay curve display Displays the full width at half maximum (FWHM), peak position and peak intensity for each profile

Multiple data loading and comparison on the same screen. Normalized processing makes multiple data comparison easy. Measurement examples Our long and proven record in fluorescence lifetime measurements is the reason our products are favored by many users in a wide range of fields.

Dye measurement of dye-sensitized solar cells





Fluorescence spectrum of chlorine compound: ex 412 nm, Φ = 0.26

The chlorine compound (1) which is a chlorophyll a (2) derivative offers great promise as a material for high-efficiency dye-sensitized solar cells. Designing molecules where electrons efficiently move from the sensitizing dye to titanium dioxide (TiO₂) requires certain optical parameters such the fluorescence lifetime and quantum yield of the dye.

Data courtesy of Prof. Dr. Hitoshi Tamiaki, Mr. Yusuke Kinoshita Laboratory of Bioorganic Chemistry Institute of Science and Engineering Ritsumeikan University

Thermally activated delayed fluorescence from organic LED material





One new technique gathering a great deal of attention for improving exciton generation efficiency in organic LED is called "Thermally Activated Delayed Fluorescence" (TADF). This is a phenomenon where delayed luminescence is observed with the usual nanosecond fluorescence components and also on the order of milliseconds from light emitted at the reverse inter system crossing from a minimum triplet exciton excitation level (T1) to a minimum singlet excitation level (S1). Making use of this thermally activated delayed fluorescence is likely to improve the external EL quantum efficiency of organic LED.

Data courtesy of Prof. Chihaya Adachi, Dr. Takane Endo, Center for Future Chemistry, Kyushu University Ayataka Endo, Mai Ogasawara, Atsushi Takahashi, Daisuke Yokoyama, Yoshimine Kato and Chihaya Adachi, *Adv. Mater.*, **21**, 4802-4906 (2009)

Measurement of fluorescence lifetime and quantum yield of stilbene compounds



Ladder type π -electron materials containing typical elements are a group of compounds with a unique electron structure reflecting the traits of the elements along with a rigid flat structure. Among these compounds, ladder type stilbene BPS possessing a phosphoryl group in cross-linked sections is unique and exhibits a light blue fluorescence in spite of a short π conjugation length as well as a high quantum yield (Φ =0.98). Moreover, measuring the fluorescence lifetime showed a result of τ =15.7 ns. From this results it is clear that, compared to stilbene compounds with cross-linked sections of other elements, these are specially unique π -electron materials possessing an extremely small nonradiative speed constant ($6.2 \times 10^7 \text{ s}^{-1}$) while maintaining a larger Stokes shift (4500 cm⁻¹).

Data courtesy of Prof. Sigehiro Yamaguchi and Dr. Aiko Fukazawa, Department of Chemistry, Graduate School of Science, Nagoya University A. Fukazawa, M. Hara, T. Okamoto, E.-C. Son, C. Xu, K. Tamao, and S. Yamaguchi, *Org. Lett.*, **10**, 913-916 (2008).

We also offer a lineup of quantum yield measurement systems <u>allowing diversified material evaluations</u> on the same sample.

Fluorescence Lifetime and Absolute PL Quantum Yield

There are two processes when substances are excited by light irradiation from the ground state to excited singlet state (S1), then deactivated to the ground state again. One is radiative process such as fluorescence and the other is a non-radiative process released as heat.

The fluorescence lifetime $\tau\,$ (tau) is defined as

$k_f + k_{nr} = 1/\tau$

where kr is the radiative rate constant and knr is the non-radiative constant.

On the other hand, the PL Quantum Yield (Φ) is expressed as the ratio of the number of photons emitted from molecules (PN_{em}) to that absorbed by molecules (PN_{abs}).

$\Phi = PN_{em} / PN_{abs}$

The PL Quantum Yield Φ is also written as

$\Phi = \mathbf{k}_{\rm f} / (\mathbf{k}_{\rm f} + \mathbf{k}_{\rm nr})$

Thus, there is a correlation between τ (tau) and Φ as shown in the following equation, and they are very important parameters for controlling the emission mechanisms of the materials.

 $k_f = \Phi / \tau$



A diversified evaluation of the luminescence materials is available!

Quantaurus-Tau for measuring fluorescence lifetime and Quantaurus-QY for absolute PL quantum yield

with simplified and minimized operating procedure are available for everybody.

Combination of Quantaurus-Tau and Quantaurus-QY allow users to obtain complementary analysis results.





Specifications

| Type number | C11367-11 | C11367-14 | C11367-12 | C11367-15 |
|-----------------------------------|--|---------------------------|---------------------|---------------------------|
| Sample | Solution, Thin-film | Solid (Thin-film, Powder) | Solution, Thin-film | Solid (Thin-film, Powder) |
| Detector type | Standard | | NIR | |
| Wavelength range | 300 nm to 800 nm | | 380 nm to 1030 nm | |
| Excitation light source | Seven types of LED light source (280 nm, 340 nm, 365 nm, 405 nm, 470 nm, 590 nm, 630 nm) | | | |
| Excitation light source switching | Software control | | | |
| Monochromator | Czerny-Turner monochromator | | | |
| Measurement time range | 2.5 ns to 50 μs / full scale | | | |
| Phosphorescence measurement | 50 μs to 50 ms / full scale | | | |
| (Option) | Phosphorescence excitation wavelength (280 nm, 340 nm, 370 nm, 405 nm, 445 nm, 470 nm, 590 nm, 630 nm) | | | |
| Time axis channel | 512 ch, 1024 ch, 2048 ch, 4096 ch | | | |
| Total time resolution | < 1.0 ns FWHM (IRF with 365 nm LED) | | | |
| Analysis function | Fluorescence lifetime analysis (up to five components by exponential function fitting) and spectrum analysis | | | |
| Supported OS | Windows 7 (32 bit) | | | |

Options

Sample box

Sample box for solution sample A12178-01

A12178-01 is a sample box for measuring the solution samples (standard: compatible with 10 mm square cells) or thin film samples. The normal sample box of C11367-11 or C11367-12 is A12178-01.

Sample box for solid sample A11551-01

A11551-01 is a sample box for measuring the powder samples or thin solid film samples. The normal sample box of C11367-14 or C11367-15 is A11551-01.

Sample box for dewar flask holder A11797-01

A11797-01 is a sample box for setting A11238-04 when measuring the lifetime of a solution samples at liquid nitrogen temperature.

Sample holder

- Dewar flask holder for low temperature measurement A11238-04
 This is used to cool the solution sample with liquid nitrogen.
- Sample holder for cryostat OPTISTAT DN This is a sample holder for setting the OPTISTAT DN (cryostat).

Sample case

For solution

- Side-arm cells (3 sets) A10095-02
- Sample tube for low temperature measurement A10095-04
 This is used to measure a sample solution at liquid nitrogen temperature.

For powder

- Laboratory dish A10095-01, -03 (with cover)
- This is used for making measurements on powder samples. This contains 5 dishes made of synthetic quartz, which suppresses fluorescence and luminescence.



Light source option

External light source connection adapter A11553-01

A11553-01 is used when excitation light source is Xe flash lamp for phosphorescence or YAG laser unit or PLP-10. It is necessary to attach A11553-01 to each of A12178-01 and A11797-01.

- External light source connection adapter A11553-02 A11553-02 is used when excitation light source is Xe flash lamp for phosphorescence or YAG laser unit or PLP-10. A11553-02 is used by attaching to A11551-01.
- Options for phosphorescence measurement C11567-01 These include xenon flash lamps, photon counting boards, and interference filter holders.
- Bandpass filter Φ25 mm for C11567-01 Selectable from among the following wavelengths: 280 nm, 340 nm, 370 nm, 405 nm, 445 nm, 470 nm, 590 nm, and 630 nm.

YAG laser unit (532 nm) C12179-01 YAG laser output power : 30 mW, Pulse width : <1.0 ns, Repetition rate : 15 kHz</p>

• YAG laser unit (355 nm) C12179-02 YAG laser output power : 4.5 mW, Pulse width : <1.0 ns, Repetition rate : 15 kHz

Picosecond Light Pulser PLP-10

The PLP-10 is picosecond light pulser which is used in combination with a controller C10196 and Laser Diode Head M12488 series. Available to select the following wavelength by replacing Laser Diode Head; 375 nm, 405 nm, 445 nm, 473 nm, 485 nm, 515 nm, 655 nm, 785 nm, 850 nm

PLP adapter A12487-01

A12487-01 is adapter for attaching M12488 to Quantaurus-Tau. A12487-01 is used when excitation light source is the PLP-10.

Dimensional outlines (unit : mm) Weight : 32 kg



Quantaurus-QY is registered trademark of Hamamatsu Photonics K.K. (China, EU, Korea, Japan, U.S.A.)

Quantaurus-Tau is registered trademark of Hamamatsu Photonics K.K. (China, EU, Korea, Japan, U.S.A.)

Product and software package names noted in this documentation are trademarks or registered trademarks of their respective manufacturers.

- Subject to local technical requirements and regulations, availability of products included in this promotional material may vary. Please consult your local sales representative.
- Information furnished by HAMAMATSU is believed to be reliable. However, no responsibility is assumed for possible inaccuracies or omissions.
- Specifications and external appearance are subject to change without notice.

© 2013 Hamamatsu Photonics K.K.

HAMAMATSU PHOTONICS K.K. www.hamamatsu.com

HAMAMATSU PHOTONICS K.K., Systems Division

812 Joko-cho, Higashi-ku, Hamamatsu City, 431-3196, Japan, Telephone: (81)53-431-0124, Fax: (81)53-435-1574, E-mail: export@sys.hpk.co.jp

U.S.A.: Hamamatsu Corporation: 360 Foothill Road, P. O. Box 6910, Bridgewater. N.J. 08807-0910, U.S.A., Telephone: (1)908-231-0960, Fax: (1)908-231-1218 E-mail: usa@hamamatsu.com

Germany: Hamamatsu Photonics Deutschland GmbH: Arzbergerstr. 10, D-82211 Herrsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-265-8 E-mail: info@hamamatsu.de France: Hamamatsu Photonics France S.A.R.L: 19, Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 10, Fax: (33)1 69 53 71 10, E-mail: info@hamamatsu.dr United Kingdom: Hamamatsu Photonics UK Limited: 2 Howard Court, 10 Tewin Road Welwyn Garden City Hertfordshire AL7 18WU, United Kingdom; Telapehone: 44(0)1707-325777 E-mail: info@hamamatsu.dr North Europe: Hamamatsu Photonics Norden AB: Torshamnsgatan 35 16440 Kista, Sweden, Telephone: (46)8-509-031-00, Fax: (46)8-509-031-01 E-mail: info@hamamatsu.se Italy: Hamamatsu Photonics (Italia: S.R.L:: Strada della Moia, 1/E, 20020 Arese, (Milano), Italy, Telephone: (39)02-935 81 733, Fax: (39)02-935 81 741 E-mail: info@hamamatsu.it China: Hamamatsu Photonics (China) Co., Ltd: 1201 Tower B, Jiaming Center, 27 Dongsanhuan Road North, Chaoyang District, Beijing 100020, China, Telephone: (86)10-6586-6006, Fax: (86)10-6586-2666 E-mail: hpc@hamamatsu.co.uk

6 E-mail: hpc@hamamatsu.com.cn Cat. No. SHSS0011E09 JUL/2013 HPK Created in Japan